Figure 1



Figure 2

## Ventilator Adaptor for use with MDI

with or without dose counter. Accepts metal or plastic tipped standard MDI canisters.

### Indications for Use:

The Southmedic MDI Adaptor is indicated for ventilated patients requiring bronchial medication.

### Intended Use:

The SMDIA-1000 MDI Adaptor is intended to be used with an MDI for intermittent delivery of prescribed aerosol medications dispensed by MDIs. It is intended to be used only when connected to ventilator tubing or tracheal tubes.

### Directions for Use:

1. Connect the SMDIA-1000 Adaptor in Line.
2. Ensure that the air flow arrows indicated on the adaptor are pointing in the direction of flow.
3. Remove the plug on the tether.
4. (a) When using an MDI with counter, refer to Figure 1.
4. (b) When using an MDI with no counter, refer to Figure 2.
5. Administer medication from the MDI canisters as prescribed.
6. When complete, replace plug.
7. Dispose of in the same manner as the circuit in which the MDI Adaptor is used.

### Warnings: ⚠

Single patient use not to be left in circuit for more than 29 days due to effect on performance, cross-contamination and risk of infection. Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury and death. Do not smoke while using this device.

### Contraindications:

Patients who cannot tolerate the addition of 16cc of dead space

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## REF SMDIA-1000

## Ventilator Adaptor for use with MDI 22mm ID x 22mm OD/15mm ID Flow ID to OD

Rx Only

Not made with natural rubber latex.

50 Alliance Blvd., Barrie, ON, Canada L4M 5K3

705-726-9383

705-728-9537 FAX

www.southmedic.com

Made in Canada

# Table 1 - Particle characterization and drug delivery results

- DE** Tabelle 1 – Partikelcharakterisierung und Medikamentenabgabeberechnung
- FR** Tableau 1 – Caractérisation des particules et résultats de la libération des médicaments
- ES** Tabla 1: Caracterización de partículas y resultados del suministro de medicamentos
- IT** Tabella 1 – Caratterizzazione delle particelle e risultati di erogazione del farmaco
- PT** Tabela 1 – Caracterização das partículas e resultados da liberação do medicamento
- SV** Tabell 1 – Partikelkaraktärisering och läkemedelsdoseringens resultat
- NL** Tabel 1 – Resultaten deeltjeskarakterisatie en medicijnafgifte
- EL** Πίνακας 1 – Χαρακτηρισμών σωματιδίων και αποτελέσματα αποδέσμευσης φαρμάκου

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A series of aerosol tests were performed on the SMDIA-1000 using an 8 stage cascade impactor at a sampling flow rate of 28 l/min equipped with a USP <601> induction port throat. The aerosol was sampled directly from the outlet of the SMDIA-1000. Table 1 is a summary of results with intervals given for a 95% confidence level.

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