Ventilator Adaptor for use with MDI

with or without dose counter. Accepts metal or plastic tipped standard MDI canisters.

Indications for Use:

The Southmedic MDI Adaptor is indicated for ventilated patients requiring bronchial medication.

Intended Use:

The SMDIA-1000 MDI Adaptor is intended to be used with an MDI for intermittent delivery of prescribed aerosol medications dispensed by MDIs. It is intended to be used only when connected to ventilator tubing or tracheal tubes.

Directions for Use:

1. Connect the SMDIA-1000 Adaptor in Line.
2. Ensure that the air flow arrows indicated on the adaptor are pointing in the direction of flow.
3. Remove the plug on the tether .
4. (a) When using an MDI with counter, refer to Figure 1.
4. (b) When using an MDI with no counter , refer to Figure 2.
5. Administer medication from the MDI canister as prescribed
6. When complete, replace plug.
7. Dispose of in the same manner as the circuit in which the MDI Adaptor is used.

Warnings: ⚠

Single patient use not to be left in circuit for more than 29 days due to effect on performance, cross-contamination and risk of infection. Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury and death. Do not smoke while using this device.

Contraindications:

Patients who cannot tolerate the addition of 16cc of dead space

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REF SMDIA-1000

Ventilator Adaptor for use with MDI 22mm ID x 22mm OD/15mm ID Flow ID to OD

Rx Only

Not made with natural rubber latex.

50 Alliance Blvd., Barrie, ON, Canada L4M 5K3

705-726-9383

705-728-9537 FAX

www.southmedic.com

Made in Canada

IFU0156L

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Table 1 - Particle characterization and drug delivery results

- DE** Tabelle 1 – Partikelcharakterisierung und Medikamentenabgabeberechnungen
- FR** Tableau 1 – Caractérisation des particules et résultats de la libération des médicaments
- ES** Tabla 1: Caracterización de partículas y resultados del suministro de medicamentos
- IT** Table 1 – Caratterizzazione delle particelle e risultati di erogazione del farmaco
- PT** Tabela 1 – Caracterização das partículas e resultados da liberação do medicamento
- SV** Tabell 1 – Partikelkaraktärisering och läkemedelsdoseringens resultat
- NL** Tabel 1 – Resultaten deeltjeskarakterisatie en medicijnafgifte
- EL** Πίνακας 1 – Χαρακτηρισμών σωματιδίων και αποτελέσματα αποδέσμευσης φαρμάκου

	Ventolin HFA	Atrovent HFA	QVAR
Total Emitted Dose per burst (ug) Summe der pro Sprühdosis abgegebenen Dosis (µg) Dose totale distribuée par unité (µg) Dosis total emitida por descarga (µg) Dose totale erogata per spruzzo (µg) Total da dosagem liberada por envio (mcg) Total Utsläpp Dosis per puff (µg) Totale afgestoten dosis per stuf (µg) Συνολική εκπνευόμενη δόση ανά ριπή (µg)	60.9 - 71.6	10.0 - 12.9	18.4 - 23.8
Particle Size (MMAD) Partikelgröße (MMAD) Dimension des particules (DAMM) Tamaño de las partículas (mediana del diametro aerodinámico de la masa) Dimensioni delle particelle (MMAD) Tamanho das Partículas (DAMM) Partikel Storlek (MMAD) Deeltjesgrootte (MMAD) Μέγεθος σωματιδίων (διάμεση αεροδυναμική διάμετρος μάζας)	1.8 - 1.9	0.9 - 1.4	0.4 - 0.6
Geometric Standard Deviation Geometrische Standardabweichung Écart type géométrique Desviación estándar geométrica Deviazione standard geometrica Desvio Padrão Geométrico Geometrisk Standard Avvikelse Geometrische standaardafwijking Γεωμετρική τυπική απόκλιση	2.4 - 2.9	4.0 - 5.1	3.1 - 4.5
Respirable Fraction (0.5-5 um) as Percent of Valve Label Inhalierbarer Anteil (0,5-5,0 µm) als Prozentangabe auf dem Ventiletikett Fraction inhalable (0.5 à 5.0 µm) exprimée comme pourcentage sur l'étiquette de la valve Fracción respirable (0.5-5.0 µm) como porcentaje de la etiqueta de la válvula Frazione respirabile (0,5-5,0 µm) in percentuale del valore di etichetta Fração Respirável (0,5-5,0 µm) Como Percentagem da marcação da válvula Andringbar Fraktion (0.5-5.0 µm) som Procent av Ventil Etikett Inandbar fraktion (0.5-5.0 µm) als percentage van ventillabel Αναπνεύσιμο κλάσμα (0,5-5 µm) ως ποσοστό που αναφέρει η ετικέτα της βαλβίδας	14.1% - 21.7%	11.1% - 15.9%	4.4% - 5.9%
Total Respirable Dose per burst (ug between 0.5-5 um) Inhalierbare Gesamtdosis pro Sprühdosis (µg zwischen 0,5 und 5,0 µm) Dose totale inhalable par unité (µg entre 0,5 et 5,0 µm) Dosis total respirable por descarga (µg entre 0,5 y 5,0 µm) Dose respirabile totale per spruzzo (µg fra 0,5 e 5,0 µm) Total da Dosagem Respirável por envio (entre 0,5 e 5,0 µm) Total Andringbar Dosis per puff (µg mellan 0,5 och 5,0 µm) Totale inandbare dosis per stuf (µg tussen 0,5 en 5,0 µm) Συνολική αναπνεύσιμη δόση ανά ριπή (µg μεταξύ 0,5 και 5 µm)	15.2 - 23.4	2.3 - 3.3	2.2 - 2.9
Medication Captured on USP Throat as Percent of Valve Label Am USP-Hals erfasste Medikation als Prozentangabe auf dem Ventiletikett Médicament capturé sur le col USP exprimé comme pourcentage sur l'étiquette de la valve Medicación capturada en la apertura USP como porcentaje de la etiqueta de la válvula Medicinalo catturato dal gomito USP in percentuale del valore di etichetta Medicamento Capturado no Canal Expressado como Percentagem na Etiqueta da Válvula Läkemedel Fångad på USP inlopp som Procent av Ventil Etikett Op USP-opening opgevangen medicijnen als percentage van ventillabel Φάρμακο που παγιδεύεται στον λαϊμό USP ως ποσοστό που αναφέρει η ετικέτα της βαλβίδας	30.8% - 42.9%	24.1% - 33.8%	24.5% - 36.8%
Medication Retained in Device as Percent of Valve Label In der Vorrichtung zurückbehaltene Medikation als Prozentangabe auf dem Ventiletikett Médicament resté dans le dispositif exprimé comme pourcentage sur l'étiquette de la valve Medicación retenida en el dispositivo como porcentaje de la etiqueta de la válvula Medicinalo rimasto nell'apparecchio in percentuale del valore di etichetta Medicamento Retido no Dispositivo Expressado como Percentagem na Etiqueta da Válvula Läkemedel kvar i Apparat som Procent av Ventil Etikett In apparat behouden medicijnen als percentage van ventillabel Φάρμακο που συγκρατείται στη συσκευή ως ποσοστό που αναφέρει η ετικέτα της βαλβίδας	30.2% - 47.7%	36.6% - 47.2%	47.3% - 58.0%
Coarse Particle Fraction (>4.7 um) as Percent of Valve Label Anteil grober Partikel (>4,7 µm) als Prozentangabe auf dem Ventiletikett Fraction particulaire grossière (>4,7 µm) exprimée comme pourcentage sur l'étiquette de la valve Fracción granular gruesa (>4,7 µm) como porcentaje de la etiqueta de la válvula Frazione di particelle grossolane (>4,7 µm) in percentuale del valore di etichetta Fração de Partícula Grossa (>4,7 µm) Expressada como Percentagem na Etiqueta da Válvula Grof Partikel Fraktion (>4,7 µm) som Procent av Ventil Etikett Grove deeltjesfractie (>4,7 µm) als percentage van ventillabel Κλάσμα αδρών σωματιδίων (>4,7 µm) ως ποσοστό που αναφέρει η ετικέτα της βαλβίδας	34.7% - 47.7%	29.4% - 48.0%	25.2% - 37.4%
Fine Particle Fraction (<4.7 um) as Percent of Valve Label Anteil feiner Partikel (<4,7 µm) als Prozentangabe auf dem Ventiletikett Fraction particulaire fine (<4,7 µm) exprimée comme pourcentage sur l'étiquette de la valve Fracción granular fina (<4,7 µm) como porcentaje de la etiqueta de la válvula Frazione di particelle fini (<4,7 µm) in percentuale del valore di etichetta Fração de Partícula Fina (<4,7 µm) Expressada como Percentagem na Etiqueta da Válvula Fin Partikel Fraktion (<4,7 µm) som Procent av Ventil Etikett Fijne-deeltjesfractie (<4,7 µm) als percentage van ventillabel Κλάσμα λεπτών σωματιδίων (<4,7 µm) ως ποσοστό που αναφέρει η ετικέτα της βαλβίδας	17.1% - 23.4%	18.0% - 23.9%	9.1% - 12.7%
Ultra-Fine Particle Fraction (<1.0 um) as Percent of Valve Label Anteil ultrafeiner Partikel (<1,0 µm) als Prozentangabe auf dem Ventiletikett Fraction particulaire ultra fine (<1,0 µm) exprimée comme pourcentage sur l'étiquette de la valve Fracción granular ultra fina (<1,0 µm) como porcentaje de la etiqueta de la válvula Frazione di particelle ultrafini (<1,0 µm) in percentuale del valore di etichetta Fração de Partícula Ultra-Fina (<1,0 µm) Expressada como Percentagem na Etiqueta da Válvula Ultra-Fin Partikel Fraktion (<1,0 µm) som Procent av Ventil Etikett Ultrafijne-deeltjesfractie (<1,0 µm) als percentage van ventillabel Κλάσμα εξαιρετικά λεπτών σωματιδίων (<1 µm) ως ποσοστό που αναφέρει η ετικέτα της βαλβίδας	6.6% - 8.0%	10.4% - 14.2%	6.3% - 9.7%

A series of aerosol tests were performed on the SMDIA-1000 using an 8 stage cascade impactor at a sampling flow rate of 28 l/min equipped with a USP <601> induction port throat. The aerosol was sampled directly from the outlet of the SMDIA-1000. Table 1 is a summary of results with intervals given for a 95% confidence level.

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