# Ventilator Adaptor for use with MDI

with or without dose counter. Accepts metal or plastic tipped standard MDI canisters.

## Indications for Use:

The Southmedic MDI Adaptor is indicated for ventilated patients requiring bronchial medication.

## Intended Use:

The SMDIA-1100 MDI Adaptor is intended to be used with an MDI for intermittent delivery of prescribed aerosol medications dispensed by MDIs. It is intended to be used only when connected to ventilator tubing or tracheal tubes.

## Directions for Use:

1. Connect the SMDIA-1100 Adaptor in Line.
2. Ensure that the air flow arrows indicated on the adaptor are pointing in the direction of flow.
3. Remove the plug on the tether .
4. (a) When using an MDI with counter, refer to Figure 1.
- (b) When using an MDI with no counter , refer to Figure 2.
5. Administer medication from the MDI canister as prescribed
6. When complete, replace plug.
7. Dispose of in the same manner as the circuit in which the MDI Adaptor is used.

## Warnings: ⚠

Single patient use not to be left in circuit for more than 29 days due to effect on performance, cross-contamination and risk of infection.

Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury and death.

Do not smoke while using this device.

## Contraindications:

Patients who cannot tolerate the addition of 16cc of dead space

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