Oxy₂Mask™

O₂ / EtCO₂

REF

Oxygen Tube Length	7 feet		14 feet	
	SLF 8"	SLM 8'	SLF 8"	SLM 15'
<i>Oxy₂Mask™ EtCO₂ Kid</i>	OK-2125-8	OK-2125-8SLM	-	-
<i>Oxy₂Mask™ EtCO₂ Adult</i>	OM-2125-8	OM-2125-8SLM	OM-2125-14	OM-2125-14SLM
<i>Oxy₂Mask™ EtCO₂ Plus</i>	OP-2125-8	OP-2125-8SLM	OP-2125-14	-

Intended Use: A minimal contact, open oxygen delivery device for use with USP medical grade oxygen for monitored breathing by providing a means to sample exhaled CO₂.

Warning: Hypercapnia will not be reliably detected by sampling through this device, the displayed CO₂ values or any other parameter/index driven from CO₂ will vary with the sampling port position, oxygen flow rate and respiratory pattern.

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Caution: Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection.

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EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Emrgo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

50 Alliance Blvd.
Barrie, Ontario
Canada L4M 5K3
705-726-9383
705-728-9537 FAX
thebetteroxygenmask.com

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

Emrgo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

Oxygen Flow Settings	Approximate O ₂ Concentrations
1 LPM	24% - 34%
2 LPM	30% - 50%
4 LPM	37% - 55%
6 LPM	47% - 70%
8 LPM	50% - 79%
10 LPM	54% - 80%
12-15 LPM	59% - 83%

*Actual FIO₂ will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect performance.
** Actual flow rate required should be prescribed in conjunction with pulse oximetry.

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