REF OKHH-1025-8 - OxyMulti-Mask™ Kid

OXYGEN MASK WITH HIGH HUMIDITY

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Intended Use: A minimal contact, open oxygen delivery device for use with USP medical grade oxygen and delivery of aerosolized drug therapy or high humidity therapy.

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Caution: Oxygen in the presence of an ignitable material and an ignition source can cause a fire resulting in severe personal injury or death. Do not smoke while using this device. Single use product not to be reused due to effect on performance, cross contamination and risk of infection.

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EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland


UK CA0086 CE 2797
Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

50 Alliance Blvd.
Barrie, Ontario
Canada L4M 5K3
705-726-9383
705-728-9537 FAX
thebetteroxygenmask.com


EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

UK CA0086 CE 2797
Emergo Consulting (UK) Limited
c/o Cr 360—UL International
Compass House, Vision Park Histon
Cambridge, CB24 9BZ
United Kingdom

*Actual FiO₂ will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect performance.
** Actual flow rate required should be prescribed in conjunction with pulse oximetry.

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Oxygen Flow Settings	Approximate O ₂ Concentrations
1 LPM	22% - 23%
2 LPM	23% - 26%
3 LPM	25% - 28%
4 LPM	28% - 33%
6 LPM	35% - 47%
8 LPM	45% - 58%
10 LPM	53% - 65%
12 LPM	58% - 71%
≥15 LPM	68% - 93%

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