



Southmedic 

Blades and Scalpels

Blades



About Southmedic

Established in 1983 Southmedic Inc. is a privately owned corporation. Southmedic specializes in product development, design, manufacturing and distribution.

Southmedic's proprietary products are sold globally through a network of distributors.

In 2008 Southmedic purchased the medical surgical blade and safety scalpel assets from Personna Medical (American Safety Razor). Southmedic updated all equipment to meet current safety and operating standards. We have successfully transitioned the branding from Personna to Southmedic Blades and Scalpels. Southmedic Blades and Scalpels are of the utmost quality.

Our constant focus in Research and Development and Innovation has allowed us to develop CABO™, the first safety scalpel cartridge to fit a standard surgical scalpel handle and complements our existing safety scalpel lines.

Polymer Coated Blades



| | | | |
|-----|--|---------|---------------------------|
| 10 | | 73-0410 | 300 Blades, 6 Boxes of 50 |
| 11 | | 73-0411 | 300 Blades, 6 Boxes of 50 |
| 12 | | 73-0312 | 50 Blades/Box |
| 15 | | 73-0415 | 300 Blades, 6 Boxes of 50 |
| 15C | | 73-4415 | 50 Blades/Box |
| 20 | | 73-0420 | 300 Blades, 6 Boxes of 50 |
| 21 | | 73-0421 | 300 Blades, 6 Boxes of 50 |
| 22 | | 73-0422 | 300 Blades, 6 Boxes of 50 |
| 23 | | 73-0323 | 50 Blades/Box |
| 24 | | 73-0424 | 50 Blades/Box |

Sterile Stainless Blades

| | | | |
|----|--|---------|----------------|
| 10 | | 73-0010 | 150 Blades/Box |
| 11 | | 73-0011 | 150 Blades/Box |
| 15 | | 73-0015 | 150 Blades/Box |
| 20 | | 73-0020 | 150 Blades/Box |
| 21 | | 73-0021 | 150 Blades/Box |
| 22 | | 73-0022 | 150 Blades/Box |

Our polymer coated blades give a smoother cut, reducing tissue drag and trauma

Our unique polymer coated edge is designed to glide with smooth, clean precision for less tissue trauma and a faster recovery. No other blade reduces drag like a Southmedic blade.



Smooth, clean cut is illustrated by cutting through 10 mm polystyrene (magnified 750 times). Our precision grinding technology coupled with our advanced polymer coating process, provides a significantly improved cutting edge.

A regular surgical blade leaves a ragged edge as indicated by a lack of uniform smoothness when cutting the polystyrene.

Non-Sterile Stainless Blades


| | | | |
|----|---|---------|-------------------|
| 22 |  | 78-0322 | 150 Blades/Carton |
| 60 |  | 78-0360 | 150 Blades/Carton |

Dermaplaning Blades

| | | | |
|--|---|-------------|---------------|
| |  | 73-0410R-50 | 50 Blades/Box |
|--|---|-------------|---------------|

Scalpels

Southmedic's disposable scalpels are available with a variety of surgical blade sizes. Designed for both right and left handed surgeons our scalpels have a depth indicator on both sides. The patented staking process of the blade provides added control, while special grooves on the handle help improve grip.

| | | | | |
|----|---|--------------------|---|--|
| 10 |  | 73-0110 73-0210 | Sterile Disposable Non-Sterile | 10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk |
| 11 |  | 73-0111 73-0211 | Sterile Disposable Non-Sterile | 10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk |
| 15 |  | 73-0115 73-0215 | Sterile Disposable Non-Sterile | 10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk |
| 20 |  | 73-0120 73-0220 | Sterile Disposable Non-Sterile | 10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk |
| 21 |  | 73-0121 73-0221 | Sterile Disposable Non-Sterile | 10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk |
| 11 |  | 73-0911 | Disposable Mini-Scalpel Non-Sterile | 5 Bags of 100 - 500/Case - Bulk |

It's the kind of difference a surgeon can actually feel.



Easy open and close

Disposable Safety Scalpels

Our disposable safety scalpels use our unique polymer coating to reduce tissue drag and trauma





- Easy open and close blade guard
- Single push lock. With just one push you can permanently lock the blade guard for safe disposal
- Large range of sterile and non-sterile sizes available
- Customized ruler in handle



| | | | | |
|----|---|--------------------|--|---|
| 10 |  | 73-0610 73-1810 | Sterile Disposable Non-Sterile Disposable | 10/Box - 12 Boxes/Case - 120/Case 50/Tray - 10 Trays/Case - 500/Case |
| 11 |  | 73-0611 73-1811 | Sterile Disposable Non-Sterile Disposable | 10/Box - 12 Boxes/Case - 120/Case 50/Tray - 10 Trays/Case - 500/Case |
| 15 |  | 73-0615 73-1815 | Sterile Disposable Non-Sterile Disposable | 10/Box - 12 Boxes/Case - 120/Case 50/Tray - 10 Trays/Case - 500/Case |



Large blade sizes for Histotechs, Pathologists, Organ Banks and Morticians

| | | | | |
|----|---|--------------------|--|---|
| 20 |  | 73-0620 73-1820 | Sterile Disposable Non-Sterile Disposable | 10/Box - 12 Boxes/Case - 120/Case 50/Tray - 10 Trays/Case - 500/Case |
| 21 |  | 73-0621 73-1821 | Sterile Disposable Non-Sterile Disposable | 10/Box - 12 Boxes/Case - 120/Case 50/Tray - 10 Trays/Case - 500/Case |
| 22 |  | 73-0622 73-1822 | Sterile Disposable Non-Sterile Disposable | 10/Box - 12 Boxes/Case - 120/Case 50/Tray - 10 Trays/Case - 500/Case |
| 60 |  | 73-1860 | Safety Autopsy Scalpel Non-Sterile Disposable | 50/Tray - 10 Trays/Case - 500/Case |

Disposable Safety Scalpels

Polymer Coated Surgical Blade

Regarded by many as the sharpest, smoothest cutting edge available

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade.

Safe Disposal

One push to lock permanently for safe disposal

Easy Open and Close

With simple locking/unlocking mechanism for safe passes

Customized Ruler

Designed into the reverse of the handle



Actual sizes shown



Larger Blade Sizes

For Histotechs, Pathologists, Organ Banks and Morticians

Disposable

Dispose of the whole scalpel safely and easily

Heavy Handle

Strong and stiff custom plastic handle. Ideal for organ retrieval and autopsy

SPSS Disposable Scalpels

Avoid sharp injuries with the first truly closed disposable safety system

- Fully retracts to expose entire blade
- SPSS cartridges are paired with Southmedic's polymer coated blade that reduces tissue drag and trauma
- Its innovative features meet OSHA Needle Stick Safety & Prevention Guidelines
- The perfect "disposable" safety scalpel system

Easy and Safe Disposal

The blade can be safely ejected from the handle reducing sharps waste and expensive disposal costs

Customized Ruler

Designed into the reverse of the handle

Polymer Coated Surgical Blade

The sharpest and smoothest cutting edge available

Safety Tab for Optimum Protection

The unique tab design effectively guards against accidental injury during installation and is completely x-ray detectable

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade

| | Sterile 10/Box 10 Boxes/Case | Non-Sterile 100/Bag 5 Bags/Case |
|-----------|---|--|
| 10 | 73-9010 | 73-9210 |
| 11 | 73-9011 | 73-9211 |
| 15 | 73-9015 | 73-9215 |
| 20 | 73-9020 | 73-9220 |
| 21 | 73-9021 | 73-9221 |
| 22 | 73-9022 | 73-9222 |



SPSS Metal Handle Safety Scalpel

The weight and balance of a regular scalpel with our unique safety system

- SPSS cartridges have Southmedic's polymer coated blade that reduces incision drag and tissue trauma
- Its innovative features meet OSHA Needle Stick Safety & Prevention Guidelines

Safety Tab for Optimum Protection

The unique tab design effectively guards against accidental injury during installation and is x-ray detectable

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade

Polymer Coated Surgical Blade

Regarded by many as the sharpest and smoothest cutting edge available

Easy Loading and Disposable Cartridge

The cartridge system provides a safe and easy way to use and dispose of blades

Reusable Handle

Cost effective reusable handle reduces waste

Customized Ruler

Designed into the reverse of the handle

Metal Handle

Handle has the same weight and feel of conventional scalpels



| | Sterile 50/Box 3 Boxes/Case | Non-Sterile 100/Bag 5 Bags/Case |
|-----------|--|--|
| 10 | 73-8010 | 73-8210 |
| 11 | 73-8011 | 73-8211 |
| 15 | 73-8015 | 73-8215 |
| 20 | 73-8020 | 73-8220 |
| 21 | 73-8021 | 73-8221 |
| 22 | 73-8022 | 73-8222 |

SPSS Safety Scalpel Metal Handles

| | | |
|-----------|-----|-------|
| 73-8030-1 | #3 | 1/pkg |
| 73-8031-1 | #3L | 1/pkg |
| 73-8040-1 | #4 | 1/pkg |
| 73-8041-1 | #4L | 1/pkg |
| 73-8071-1 | #7 | 1/pkg |



CABO

Cutting Above and Beyond Objection

The world's only safety scalpel that fits the surgeon's favorite handle

- Its innovative features meet OSHA Needle Stick Safety & Prevention Guidelines



Easy Open & Close

With simple locking/unlocking mechanism

Easy Loading and Disposable

The cartridge system provides a safe and easy way to use and dispose of blades and is completely x-ray detectable

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade

Polymer Coated Surgical Blade

The sharpest and smoothest cutting edge available - reduces tissue drag and trauma

Reusable Handle

The most cost effective Safety Scalpel available on the market. Reusable handles create less disposable waste



| | | | | | | |
|----|---|---------|------------------------|---------|--------------|----------|
| 10 |  | 73-7010 | Sterile Disposable | 50/Box | 3 Boxes/Case | 150/Case |
| | | 73-7210 | Non-Sterile Disposable | 100/Box | 5 Boxes/Case | 500/Case |
| 11 |  | 73-7011 | Sterile Disposable | 50/Box | 3 Boxes/Case | 150/Case |
| | | 73-7211 | Non-Sterile Disposable | 100/Box | 5 Boxes/Case | 500/Case |
| 15 |  | 73-7015 | Sterile Disposable | 50/Box | 3 Boxes/Case | 150/Case |
| | | 73-7215 | Non-Sterile Disposable | 100/Box | 5 Boxes/Case | 500/Case |
| 20 |  | 73-7020 | Sterile Disposable | 50/Box | 3 Boxes/Case | 150/Case |
| | | 73-7220 | Non-Sterile Disposable | 100/Box | 5 Boxes/Case | 500/Case |
| 21 |  | 73-7021 | Sterile Disposable | 50/Box | 3 Boxes/Case | 150/Case |
| | | 73-7221 | Non-Sterile Disposable | 100/Box | 5 Boxes/Case | 500/Case |
| 22 |  | 73-7022 | Sterile Disposable | 50/Box | 3 Boxes/Case | 150/Case |
| | | 73-7222 | Non-Sterile Disposable | 100/Box | 5 Boxes/Case | 500/Case |



Manufactured by:
Southmedic Inc.
50 Alliance Blvd., Barrie, ON L4M 5K3
www.southmedic.com
smedic@southmedic.com

1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
ISO 13485