

# Endotracheal Tube Protector

## with Security Ring



The Southmedic Bite Block protects endotracheal tubes (ET tubes) when sedated patients bite down

- Prevents occlusion of the ET tube
- Fits multiple ET tube sizes
- Soft enough to protect teeth
- Transparent material aids tube position
- Security rings fit over standard ET tubes to prevent slipping
- Tie-wrap to secure Bite Block
- Does not contain natural rubber latex or phthalates such as DEHP
- Available in two standard sizes
- Ready to use
- For single patient use
- 50 per case

### New design:

- Triple security rings increase tether length, allowing greater bite block placement control
- New one piece design, incorporates built-in adjustment strap; making it easier to attach and adjust
- Improved design accommodates most ET tube attachment devices and fasteners



**New** ETPSR3-00A50  
Adult: 6.5mm-9.5mm  
Triple Security Ring



ETPSR-00A50  
Adult: 6.5mm-9.5mm  
Single Security Ring



ETPSR-00P50  
Pediatric 3.5mm-6.0mm  
Single Security Ring



Manufactured by:  
Southmedic Inc.  
50 Alliance Blvd., Barrie, ON L4M 5K3  
www.southmedic.com  
contactus@southmedic.com

1-800-463-7146 Toll free in North America  
1-705-726-9383  
1-705-728-9537 Fax  
ISO 13485