

MDI Adaptor for Dosage Counters

Southmedic inline MDI Adaptors ensure safe and complete delivery of prescribed aerosol medication

- Compatible with plastic and metal canister stems
- Can be used with, or without dose counter
- Ready to use
- For single patient use
- Includes built in tether strap
- Does not contain natural latex and DEHP/Phthalates
- Available in four standard sizes
- 25 per case

Southmedic SMDI adaptors remain inline:

- Improves patient safety; reduced opportunity for cross-contamination
- Improves patient comfort; reduced manipulation of circuit
- Improves compliance of prescribed medication due to reduced interruptions of care



SMDIA-1000



SMDIA-1100

NOTE: Arrows indicate direction of flow



SMDIA-1000
22mmID → 22mmOD/15mmID



SMDIA-1100
22mmOD → 22mmID



SMDIA-1200
22mmOD → 22mmOD/15mmID



SMDIA-1300
15mmOD → 15mmID



Manufactured by:
Southmedic Inc.
50 Alliance Blvd., Barrie, ON L4M 5K3
www.southmedic.com
contactus@southmedic.com

1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
ISO 13485