

CO₂ Sampling and Oxygen Delivery Cannula



- Side Stream CO₂ Sampling with simultaneous oxygen delivery
- Soft anatomically curved nasal prongs for a comfortable fit
- Anatomically designed block provides a stable fit without twisting

CO₂ is collected at the nares with a supplemental port to collect orally, while supplemental O₂ is delivered upwards to the nostrils and downwards in the direction of the mouth.

25 units/box Description

ADULT

NC-2025-6	CO ₂ /O ₂ nasal cannula with O ₂ barb connector and CO ₂ connector female luer
NC-2025-7	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and CO ₂ connector female luer
NC-2025-8F	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line female luer
NC-2025-8	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line male luer
NC-2025-14	CO ₂ /O ₂ nasal cannula with 14' O ₂ with universal connector and 14' CO ₂ sampling line male luer

Paediatric

NCP-2025-6	CO ₂ /O ₂ nasal cannula with O ₂ barb connector and CO ₂ connector female luer
NCP-2025-7	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and CO ₂ connector female luer
NCP-2025-8	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line male luer
NCP-2025-14	CO ₂ /O ₂ nasal cannula with 14' O ₂ with universal connector and 14' CO ₂ sampling line male luer



705-726-9383 smedic@southmedic.com
50 Alliance Blvd., Barrie, ON., Canada L4M 5K3

www.southmedic.com
ISO 13485

www.oxyarm.com