

Vapofil™

Vaporizer Filler Tubes That Work First Drop to Last



With Stronger and More Reliable Features


Features:

- Flexible outer tubing
- Separate inner tubes for air and agent
- Inner sleeve reinforced by brass tubing
- Colour coded and keyed for agent specific use
- Anti-pollution plug
- Durable
- Rapid filling

Ordering Info

Vaporizer	Cat. No.
■ Florane/Isoflurane	8907-F
■ Halothane/Fluothane	8907-H
■ Ultane/Sevoflurane	8907-S

Vapofil™



AntiSpil™



Features:

- Curved tip and flow meter let you pour the agent with less risk of spills
- Colour coded and keyed cap for agent specific use
- Anti-pollution plug
- Durable

Ordering Info

Vaporizer	Cat. No.
■ Florane/Isoflurane	9010-F
■ Halothane/Fluothane	9010-H
■ Ultane/Sevoflurane	9010-S



705-726-9383

www.southmedic.com
ISO 13485

email: smedic@southmedic.com

50 Alliance Blvd., Barrie, ON., Canada L4M 5K3