

Foam Head Support

Elevates and supports patient's head and neck



Angled Head Support



Circular Head Support

- Centre aperture allows access for endotracheal tube when patient is in the prone position
- Soft foam
- Single patient use

Code	Description	Centre Aperture	Dimensions (LxWxH)	Case Qty
FMH-1000-L	Angled Head Support Large	3.88"	9.00" x 7.75" x 3.85"	1
FMH-1000-S	Angled Head Support Small	3.25"	7.50" x 6.50" x 1.60"	1
FMD-1006-L	Circular Head Support Large	4"	9.00" x 9.00" x 1.50"	6
FMD-1006-S	Circular Head Support Small	2.50"	6.00" x 6.00" x 2.50"	6



Manufactured by:
Southmedic Inc.
50 Alliance Blvd., Barrie, ON L4M 5K3
www.southmedic.com
contactus@southmedic.com

1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
ISO 13485