

Blades & Scalpels



Southmedic 










Blades









About Southmedic

Established in 1983 Southmedic Inc. is a privately owned corporation. Southmedic specializes in product development, design, manufacturing and distribution. Southmedic's proprietary products are sold globally through a network of distributors.

Polymer Coated Blades



10		73-0410	300 Blades, 6 Boxes of 50
11		73-0411	300 Blades, 6 Boxes of 50
12		73-0312	50 Blades/Box
15		73-0415	300 Blades, 6 Boxes of 50
15C		73-4415	50 Blades/Box
20		73-0420	300 Blades, 6 Boxes of 50
21		73-0421	300 Blades, 6 Boxes of 50
22		73-0422	300 Blades, 6 Boxes of 50
23		73-0323	50 Blades/Box

Sterile Stainless Blades

10		73-0010	150 Blades/Box
11		73-0011	150 Blades/Box
15		73-0015	150 Blades/Box
20		73-0020	150 Blades/Box
21		73-0021	150 Blades/Box
22		73-0022	150 Blades/Box

Our polymer coated blades give a smoother cut, reducing tissue drag and trauma

Our unique polymer coated edge is designed to glide with smooth, clean precision for less tissue trauma and a faster recovery. No other blade reduces drag like a Southmedic blade.



Smooth, clean cut is illustrated by cutting through 10 mm polystyrene (magnified 750 times). Our precision grinding technology coupled with our advanced polymer coating process, provides a significantly improved cutting edge.

A regular surgical blade leaves a ragged edge as indicated by a lack of uniform smoothness when cutting the polystyrene.

Non-Sterile Stainless Blades

22		78-0322	150 Blades/Carton
60		78-0360	150 Blades/Carton

Dermaplaning Blades

	73-0410R-50	50 Blades/Box
---	-------------	---------------

Scalpels

Southmedic's disposable scalpels are available with a variety of surgical blade sizes. Designed for both right and left handed surgeons our scalpels have a depth indicator on both sides. The patented staking process of the blade provides added control, while special grooves on the handle help improve grip.

10		73-0110 73-0210	Sterile Disposable Non-Sterile	10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk
11		73-0111 73-0211	Sterile Disposable Non-Sterile	10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk
15		73-0115 73-0215	Sterile Disposable Non-Sterile	10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk
20		73-0120	Sterile Disposable	10 Scalpels/Box - 10 Boxes/Case
21		73-0121 73-0221	Sterile Disposable Non-Sterile	10 Scalpels/Box - 10 Boxes/Case 1,000 Scalpels/Case - Bulk
11		73-0911	Disposable Mini-Scalpel Non-Sterile	5 Bags of 100 - 500/Case - Bulk

It's the kind of difference a surgeon can actually feel.



Easy open and close

Disposable Safety

Our disposable safety scalpels use our unique polymer coating to reduce tissue drag and trauma

- Easy open and close blade guard
- Single push lock. With just one push you can permanently lock the blade guard for safe disposal
- Large range of sterile and non-sterile sizes available
- Customized ruler in handle



10		73-0610	Sterile Disposable	120/Case - 10/Box - 12 Boxes/Case (12 Boxes of 10)
		73-1810	Non-Sterile Disposable	500/Case - 50/Tray - 10 Trays/Case (10 Trays of 50)
11		73-0611	Sterile Disposable	120/Case - 10/Box - 12 Boxes/Case (12 Boxes of 10)
		73-1811	Non-Sterile Disposable	500/Case - 50/Tray - 10 Trays/Case (10 Trays of 50)
15		73-0615	Sterile Disposable	120/Case - 10/Box - 12 Boxes/Case (12 Boxes of 10)
		73-1815	Non-Sterile Disposable	500/Case - 50/Tray - 10 Trays/Case (10 Trays of 50)



Large blade sizes for Histotechs, Pathologists, Organ Banks and Morticians

20		73-0620	Sterile Disposable	120/Case - 10/Box - 12 Boxes/Case (12 Boxes of 10)
		73-1820	Non-Sterile Disposable	500/Case - 50/Tray - 10 Trays/Case (10 Trays of 50)
21		73-0621	Sterile Disposable	120/Case - 10/Box - 12 Boxes/Case (12 Boxes of 10)
		73-1821	Non-Sterile Disposable	500/Case - 50/Tray - 10 Trays/Case (10 Trays of 50)
22		73-0622	Sterile Disposable	120/Case - 10/Box - 12 Boxes/Case (12 Boxes of 10)
		73-1822	Non-Sterile Disposable	500/Case - 50/Tray - 10 Trays/Case (10 Trays of 50)
60		73-1860	Safety Autopsy Scalpel Non-Sterile Disposable	500/Case - 50/Tray - 10 Trays/Case (10 Trays of 50)

Scalpels

Polymer Coated Surgical Blade

Regarded by many as the sharpest, smoothest cutting edge available

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade.

Safe Disposal

One push to lock permanently for safe disposal

Easy Open and Close

With simple locking/unlocking mechanism for safe passes

Customized Ruler

Designed into the reverse of the handle



Actual sizes shown



Larger Blade Sizes

For Histotechs, Pathologists, Organ Banks and Morticians

Disposable

Dispose of the whole scalpel safely and easily

Heavy Handle

Strong and stiff custom plastic handle. Ideal for organ retrieval and autopsy

SPSS Disposable Scalpels

Avoid sharps injuries with the first truly closed disposable safety system

- Fully retracts to expose entire blade
- SPSS cartridges are paired with Southmedic's polymer coated blade that reduces tissue drag and trauma
- Its innovative features meet OSHA Needle Stick Safety & Prevention Guidelines
- The perfect "disposable" safety scalpel system

Easy and Safe Disposal

The blade can be safely ejected from the handle reducing sharps waste and expensive disposal costs

Customized Ruler

Designed into the reverse of the handle

Polymer Coated Surgical Blade

The sharpest and smoothest cutting edge available


Safety Tab for Optimum Protection

The unique tab design effectively guards against accidental injury during installation and is completely x-ray detectable

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade

	Sterile 10/Box 10 Boxes/Case	Non-Sterile 100/Bag 5 Bags/Case
10	73-9010	73-9210
11	73-9011	73-9211
15	73-9015	73-9215
20	73-9020	73-9220
21	73-9021	73-9221
22	73-9022	73-9222



SPSS Metal Handle Safety Scalpel

The weight and balance of a regular scalpel with our unique safety system

- SPSS cartridges have Southmedic's polymer coated blade that reduces tissue drag and tissue trauma
- Its innovative features meet OSHA Needle Stick Safety & Prevention Guidelines

Safety Tab for Optimum Protection

The unique tab design effectively guards against accidental injury during installation and is x-ray detectable

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade

Polymer Coated Surgical Blade

Regarded by many as the sharpest and smoothest cutting edge available

Easy Loading and Disposable Cartridge

The cartridge system provides a safe and easy way to use and dispose of blades

Reusable Handle

Cost effective reusable handle reduces waste

Customized Ruler

Designed into the reverse of the handle

Metal Handle

Handle has the same weight and feel of conventional scalpels



	Sterile 50/Box 3 Boxes/Case	Non-Sterile 100/Bag 5 Bags/Case
--	--	--

10	73-8010	73-8210
11	73-8011	73-8211
15	73-8015	73-8215
20	73-8020	73-8220
21	73-8021	73-8221
22	73-8022	73-8222

SPSS Safety Scalpel Metal Handles

73-8030-1	#3	1/pkg
73-8031-1	#3L	1/pkg
73-8041-1	#4L	1/pkg
73-8071-1	#7	1/pkg



CABO

Cutting Above and Beyond Objection

The world's only safety scalpel that fits the surgeon's favorite handle

- Its innovative features meet OSHA Needle Stick Safety & Prevention Guidelines



Easy Open & Close

With simple locking/unlocking mechanism

Easy Loading and Disposable

The cartridge system provides a safe and easy way to use and dispose of blades and is completely x-ray detectable

Clear Plastic Sheath

The design offers complete visibility and incorporates a convenient thumb rest for ease of use. Retracts fully to expose entire blade

Polymer Coated Surgical Blade

The sharpest and smoothest cutting edge available - reduces tissue drag and trauma

Reusable Handle

The most cost effective Safety Scalpel available on the market. Reusable handles create less disposable waste



10		73-7010	Sterile Disposable	50/Box	3 Boxes/Case	150/Case
		73-7210	Non-Sterile Disposable	100/Box	5 Boxes/Case	500/Case
11		73-7011	Sterile Disposable	50/Box	3 Boxes/Case	150/Case
		73-7211	Non-Sterile Disposable	100/Box	5 Boxes/Case	500/Case
15		73-7015	Sterile Disposable	50/Box	3 Boxes/Case	150/Case
		73-7215	Non-Sterile Disposable	100/Box	5 Boxes/Case	500/Case
20		73-7220	Non-Sterile Disposable	100/Box	5 Boxes/Case	500/Case
21		73-7021	Sterile Disposable	50/Box	3 Boxes/Case	150/Case
		73-7221	Non-Sterile Disposable	100/Box	5 Boxes/Case	500/Case
22		73-7022	Sterile Disposable	50/Box	3 Boxes/Case	150/Case
		73-7222	Non-Sterile Disposable	100/Box	5 Boxes/Case	500/Case



Blades & Scalpels

Respiratory

Surgical

Dynamic Tissue Systems

Manufactured by:
Southmedic Inc.
50 Alliance Blvd.,
Barrie, ON
Canada L4M 5K3
ISO 13485

1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
contactus@southmedic.com
www.southmedic.com

P0079 Mar 2018