

# MDI Adaptor for Dosage Counters

Southmedic inline MDI Adaptors ensure safe and complete delivery of prescribed aerosol medication

- Compatible with plastic and metal canister stems
- Can be used with, or without dose counter
- Ready to use
- For single patient use
- Includes built in tether strap
- Not made with natural latex and DEHP/Phthalates
- Available in four standard sizes
- 25 per case

Southmedic MDI adaptors remain inline:

- Improves patient safety; reduced opportunity for cross-contamination
- Improves patient comfort; reduced manipulation of circuit
- Improves compliance of prescribed medication due to reduced interruptions of care



SMDIA-1000



SMDIA-1100



SMDIA-1000  
22mmID -> 22mmOD/15mmID



SMDIA-1100  
22mmOD -> 22mmID



SMDIA-1200  
22mmOD -> 22mmOD/15mmID



SMDIA-1300  
15mmOD -> 15mmID

NOTE: Arrows indicate direction of flow



Respiratory

Surgical

Blades & Scalpels

Dynamic Tissue Systems

Manufactured by:  
Southmedic Inc.  
50 Alliance Blvd.,  
Barrie, ON Canada L4M 5K3  
www.southmedic.com

email: [contactus@southmedic.com](mailto:contactus@southmedic.com)  
1-800-463-7146 Toll free in North America  
1-705-726-9383  
1-705-728-9537 Fax  
ISO 13485

P0169 Feb 2018