

CO₂ Sampling and Oxygen Delivery Cannula



- Side Stream CO₂ Sampling with simultaneous oxygen delivery
- Soft anatomically curved nasal prongs for a comfortable fit
- Anatomically designed block provides a stable fit without twisting
- CO₂ is collected at the nares with a supplemental port to collect orally, while supplemental O₂ is delivered upwards to the nostrils and downwards in the direction of the mouth
- 25 per box

| Code | Description | Size | Tubing Length | Luer Type |
|-------------|--|------------|---------------|-----------|
| NC-2025-7 | CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and CO ₂ connector female luer | Adult | 7' | Female |
| NC-2025-8F | CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line female luer | Adult | 7' | Female |
| NC-2025-8 | CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line male luer | Adult | 7' | Male |
| NCD-2025-8* | CO ₂ /O ₂ divided cannula with 7' O ₂ with CO ₂ lines | Adult | 7' | Male |
| NC-2025-14 | CO ₂ /O ₂ nasal cannula with 14' O ₂ with universal connector and 14' CO ₂ sampling line male luer | Adult | 14' | Male |
| NCP-2025-8 | CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line male luer | Paediatric | 7' | Male |

*Available in Canada Only

Southmedic

Respiratory

Surgical

Blades & Scalpels

Dynamic Tissue Systems

Manufactured by:
Southmedic Inc.
50 Alliance Blvd.,
Barrie, ON
Canada L4M 5K3
ISO 13485

www.southmedic.com
email: contactus@southmedic.com
1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax

P0076 Mar 2018