

# Vapofil™


Vaporizer Filler Tubes That Work First Drop to Last



With stronger and more reliable features

- Flexible outer tubing
- Separate inner tubes for air and agent
- Inner sleeve reinforced by brass tubing
- Colour coded and keyed for agent specific use
- Anti-pollution plug
- Durable
- Rapid filling

## Vapofil™



## AntiSpil™



- Curved tip and flow meter let you pour the agent with less risk of spills
- Colour coded and keyed cap for agent specific use
- Anti-pollution plug
- Durable

### Ordering Info

Vaporizer	Cat. No.
Florane/Isoflurane	8907-F
Halothane/Fluothane	8907-H
Ultane/Sevoflurane	8907-S

### Ordering Info

Vaporizer	Cat. No.
Florane/Isoflurane	9010-F
Halothane/Fluothane	9010-H
Ultane/Sevoflurane	9010-S

**Southmedic** 

Surgical  
Respiratory  
Blades & Scalpels  
Dynamic Tissue Systems

Manufactured by:  
Southmedic Inc.  
50 Alliance Blvd.,  
Barrie, ON  
Canada L4M 5K3  
ISO 13485

www.southmedic.com  
email: [contactus@southmedic.com](mailto:contactus@southmedic.com)  
1-800-463-7146 Toll free in North America  
1-705-726-9383  
1-705-728-9537 Fax

P0059 Apr 2018