

Southmedic

DETERGENTS



- Non-Toxic
- 100% Biodegradable
- No Dyes
- 100% Water Soluble
- Non-Corrosive
- Low-Foaming
- Rust Inhibitors
- Neutral pH



2.5 Gallon
Container Shown

EcoZyme

ULTRA

MULTI-TIERED ENZYMATIC DETERGENT / DÉTERGENT MULTI-ENZYMATIQUE

EcoZyme Ultra is an extremely effective multi-enzymatic detergent for cleaning metal, plastic, dental, and veterinarian surgical instrumentation including flexible & rigid endoscopes. EcoZyme's unique patented formulation contains a variety of safe enzymes and surfactants that break down organic and inorganic biological residue.

Part ID	Volume	UOM
18560-04	1 gal	4/Box
18560-2.5	2.5 gal	2/Box
18560-05	5 gal	Each
18560-15	15 gal	Each

EcoBase

NON-ENZYMATIC DETERGENT / NON-DÉTERGENT MULTI-ENZYMATIQUE

EcoBase Detergent is a non-enzymatic detergent designed for automated processes, such as tunnel washer, ultrasonic washers and cart washers. It is an effective detergent used for the cleaning of metallic, plastic, surgical and dental instruments, including flexible and rigid endoscopes.

Part ID	Volume	UOM
11050-04	1 gal	4/Box
11050-2.5	2.5 gal	2/Box
11050-05	5 gal	Each
11050-15	15 gal	Each



15 Gallon
Container Shown



EcoZyme

SPRAY

Our same patented formula as EcoZyme Ultra but in a spray bottle is available in 8, 16, and 32 fluid ounces. EcoZyme is safe on all surgical instruments; flexible and rigid scopes and power tools. EcoZyme's unique patented formulation contains a variety of safe enzymes and surfactants designed to break down all organic and inorganic residue.



32 fl. oz. Spray Bottle Shown

Part ID	Volume	UOM
27545-08	8 oz.	12/Box
27545-16	16 oz.	12/Box
27545-32	32 oz.	6/Box



1 Gallon Container Shown

EcoZyme

AUTO

MULTI-TIERED ENZYMATIC DETERGENT / DÉTERGENT MULTI-ENZYMATIQUE

EcoZyme Auto is an effective multi-enzymatic detergent for cleaning metal, plastic, dental, veterinarian and surgical instrumentation. EcoZyme's unique patented formulation contains a variety of safe enzymes and surfactants that break down the compounds of organic and inorganic biological residue.

Part ID	Volume	UOM
17560-04	1 gal	4/Box
17560-2.5	2.5 gal	2/Box
17560-05	5 gal	Each

EcoSonic

DESCALER

EcoSonic Descaler detergent is an advanced surfactant cleaner formulated to refurbish surgical instrumentation without the use of hazardous ingredients. EcoSonic is excellent for removing stains, discolouration, rust, and residue from metal and porous/metallic surfaces without negative impact. EcoSonic prevents redistribution of contaminants. The corrosion inhibitor leaves surfaces clean, bright and residue free.

Part ID	Volume	UOM
43394-04	1 gal	4/Box



1 Gallon Container Shown



Finding a Better Way

Scan the QR code to learn more about Southmed's detergent products online.

