

Clear Choice™

Face and Eye Splash Protection

- Cost effective, re-usable plastic frames
- Does not contain BPA or natural rubber latex
- Comfortably worn over glasses and surgical masks
- OSHA Bloodborne Pathogen Standard and EEC Mandate Compliant 1910.1030(d)93(x)



	No. of Frames Included	No. of Shields Included	Shield Size in inches (L x H x D)	Securing Mechanism	Case Pack
--	------------------------	-------------------------	-----------------------------------	--------------------	-----------


Clear Choice™ Face Shields

SMI45100S-PF	-	100	9.75 x 8.56 x .0065	-	100
SMI2510-PF	-	25	9.75 x 8.56 x .0065	-	25
SMI1010	10	10	9.75 x 8.56 x .0065	Adjustable Elastic Tube	10



Clear Choice™ Eye Shields

SMI9210-100	-	100	9.0 x 2.625 x .006	-	100
SMI9212	100	100	9.0 x 2.625 x .006	Non-Folding Arms	100
SMI9452	10	25	9.0 x 2.625 x .006	Non-Folding Arms	10



Clear Choice™ Plastic Frames

SMI9211MC	100	-	-	Non-Folding Arms	100
SMI9214MC	100	-	-	Folding Arms	100

Clear Choice™ Eye Shield Dispenser Box

SMI9001	-	-	-	-	1
---------	---	---	---	---	---



Clear Choice™ Foam Face Shield

comfortable foam head rest • two shield lengths • single use • non-adjustable elastic

SMI3610	-	40	13.0 x 7.5 x .006	Elastic Headband	40
SMI3611	-	40	13.0 x 9.5 x .006	Elastic Headband	40

Southmedic

Surgical

Respiratory

Blades & Scalpels

Dynamic Tissue Systems

Manufactured by:
Southmedic Inc.
50 Alliance Blvd.,
Barrie, ON
Canada L4M 5K3
ISO 13485

1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
contactus@southmedic.com
www.southmedic.com

Nov 2017 P001107