


# Foam Head Support

Elevates and supports patient's head and neck



Angled Head Support



Circular Head Support

- Centre aperture allows access for endotracheal tube when patient is in the prone position
- Soft foam
- Single use

Code	Description	Centre Aperture	Dimensions (LxWxH)	Case Qty
FMH-1000-L	Angled Head Support Large	3.88"	9.00" x 7.75" x 3.85"	1
FMH-1000-S	Angled Head Support Small	3.25"	7.50" x 6.50" x 1.60"	1
FMD-1006-L	Circular Head Support Large	4"	9.00" x 9.00" x 1.50"	6



Surgical  
Respiratory  
Blades & Scalpels  
Dynamic Tissue Systems

Manufactured by:  
Southmedic Inc.  
50 Alliance Blvd.,  
Barrie, ON  
Canada L4M 5K3  
ISO 13485

1-800-463-7146 Toll free in North America  
1-705-726-9383  
1-705-728-9537 Fax  
contactus@southmedic.com  
www.southmedic.com

Nov/2017 P0209