

# OxyMask™

Oxygen Therapy



## Patient Safety

- Reduces the risk of CO<sub>2</sub> rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather
- Available EtCO<sub>2</sub> monitoring

## Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

## Oxygen Therapy Compliance


Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy

## Operational Efficiencies

- One OxyMask™ replaces multiple traditional oxygen therapy interfaces:
- SKU standardization
  - Inventory reduction
  - Reduced medical waste
  - Reduces complexity for bedside staff

# OxyMask™

- Delivers Oxygen flows from 1 to 15+ litres per minute - flush, 24% to 90% FiO<sub>2</sub>  
(Tyke size 1/4L to 5+, 22% to 65% FiO<sub>2</sub>)
- Four sizes available:
  - Adult
  - Plus (25% larger mask)
  - Kid (3 - 10 years old, 33 - 70 lbs / 15 - 32 kg)
  - Tyke (6 months - 3 years old, 16 - 33 lbs / 7 - 15 kg)
- Does not contain natural rubber latex or phthalates such as DEHP
- 25 masks/case



Code	Description	Mask Size	Tubing Length	Elastic
OM-1125-8	OxyMask™ with 7' tubing	Adult	7'	Single
OM-1125-14	OxyMask™ with 14' tubing	Adult	14'	Single
OM-1425-8	OxyMask™ with ear elastics and 7' tubing	Adult	7'	Double
OM-7025-8*	OxyMask™ Non PVC with tubing 7'	Adult	7'	Single
OP-1125-8	OxyMask™ Plus w 7' tubing	Plus	7'	Single
OP-1125-14	OxyMask™ Plus w 14' tubing	Plus	14'	Single
OP-1425-8	OxyMask™ Plus w ear elastics & 7' tubing	Plus	7'	Double
OP-7025-8*	OxyMask™ Plus Non PVC with tubing 7'	Plus	7'	Single
OK-1125-8	OxyMask™ Kid w 7' tubing	Kid	7'	Single
OK-1125-14	OxyMask™ Kid w 14' tubing	Kid	14'	Single
OK-7025-8*	OxyMask™ Kid Non PVC with tubing 7'	Kid	7'	Single
OT-1025-8	OxyMask™ Tyke w 7' tubing	Tyke	7'	Single
OT-7025-8*	OxyMask™ Tyke Non PVC with tubing 7'	Tyke	7'	Single

\*Non PVC



Manufactured by:  
Southmedic Inc.  
50 Alliance Blvd.,  
Barrie, ON Canada L4M 5K3  
www.thebetteroxygenmask.com

email: [oxymask@southmedic.com](mailto:oxymask@southmedic.com)  
1-800-463-7146 Toll free in North America  
1-705-726-9383  
1-705-728-9537 Fax  
ISO 13485

Jan 2018 P0177