


Foam Head Support

Elevates and supports patient's head and neck



- Centre aperture allows access for endotracheal tube when patient is in the prone position
- Soft foam
- Single use

Code	Description	Centre Aperture	Dimensions (LxWxH)	Case Qty
FMH-1000-L	Angled Head Support Large	3.88"	9.00" x 7.75" x 3.85"	1
FMH-1000-S	Angled Head Support Small	3.25"	7.50" x 6.50" x 1.60"	1



Surgical
Respiratory
Blades & Scalpels
Dynamic Tissue Systems

Manufactured by:
Southmedic Inc.
50 Alliance Blvd.,
Barrie, ON
Canada L4M 5K3
ISO 13485

1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
contactus@southmedic.com
www.southmedic.com

Mar 2018 P0209