

Endotracheal Tube Protector

with Security Ring



The Southmedic Bite Block protects endotracheal tubes (ET tubes) when sedated patients bite down

- Prevents occlusion of the ET tube
- Fits multiple ET tube sizes
- Soft enough to protect teeth
- Transparent material aids tube position
- Security rings fit over standard ET tubes to prevent slipping
- Tie-wrap to secure Bite Block
- Not made with natural rubber latex or phthalates such as DEHP
- Available in two standard sizes
- Ready to use
- For single use
- 50 per case

New design:

- Triple security rings increase tether length, allowing greater bite block placement control
- New one piece design, incorporates built-in adjustment strap; making it easier to attach and adjust
- Improved design accommodates most ET tube attachment devices and fasteners



New ETPSR3-00A50
Adult: 6.5mm-9.5mm
Triple Security Ring



ETPSR-00A50
Adult: 6.5mm-9.5mm
Single Security Ring



ETPSR-00P50
Pediatric 3.5mm-6.0mm
Single Security Ring

Southmedic

Respiratory

Surgical

Blades & Scalpels

Dynamic Tissue Systems

Manufactured by:
Southmedic Inc.
50 Alliance Blvd.,
Barrie, ON Canada L4M 5K3
www.southmedic.com

email: contactus@southmedic.com
1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax
ISO 13485

P0170 Feb 2018