

EcoZyme

SPRAY



Southmedic’s EcoZyme Spray; the only brand that can promise “*Every Ounce is Safe*”! Our spray is designed to break down all bioburden, body fluids and blood on contact, preventing it from drying on instruments before they make their way down to the MDR department. Our spray is also ideal for flexible or rigid scopes and power tools.



As stated by various governing bodies such as CSA, ORNAC or AAMI: saline solutions, abrasive chemicals, or harsh cleaners shall not be used. We’ve got you covered! Our EcoZyme Spray has ZERO hazardous chemicals, is odourless, 100% biodegradable, non-corrosive, and non-toxic; making it safe for it’s users and all instrumentation it is applied to.

Our purpose is to provide the safest cleaning detergents for hospital staff, to minimize or eliminate the problems that arise while using hazardous materials, creating a more pleasant work environment.

PART ID: #27545-04

4 oz.
72 Units / Box
**not shown*

PART ID: #27545-08

8 oz.
12 Units / Box

PART ID: #27545-16



16 oz.
12 Units / Box

PART ID: #27545-32

32 oz.
6 Units / Box

We recommend wearing PPE, keeping instruments in the open position, spraying 6-8” from the heavily soiled instruments once in their container or case. Making sure all heavily contaminated instruments are completely coated and close or cover the container prior to transportation. The Ecozyme spray starts working right away!

EcoZyme Spray Designed for: Point-of-Use Pre-Spray



www.southmedic.com