

EcoZyme

ULTRA



Southmedic's EcoZyme Ultra; the only brand that can promise “*Every Ounce is Safe*”! EcoZyme Ultra is an extremely effective multi-enzymatic detergent for cleaning metal, plastic, dental, and veterinarian surgical instrumentation including flexible & rigid endoscopes.



- Highly concentrated formula means you don't need as much product
- EcoZyme starts lifting debris and soils right away!
- Less scrubbing and manual cleaning; allowing technicians to reprocess instrumentation quickly
- Our clear amber colour, and low foam in the sink allows for clarity and less opportunity for injury
- The non corrosive, and rust inhibitor properties of EcoZyme does not harm instruments and devices; protecting them from damage, even after repeated or extended exposure; keeping cost to replace instruments down
- Stable supply chain! Our Detergents are manufactured, sourced, and warehoused locally here in Ontario, Canada, and stocked in Illinois to service all of North America efficiently
- No contributing VOC's (volatile organic compounds)
- Our detergents are odourless and 100% biodegradable

PART ID: #18560-04

1 Gallon
4 Units / Box

PART ID: #18560-2.5

2.5 Gallons
2 Units / Box

PART ID: #18560-05

5 Gallons
Each

PART ID: #18560-15


15 Gallons
Each
**not shown*

Our EcoZyme Ultra detergent contains specific enzymes that quickly and efficiently dissolve soils off instruments. It also contains sequestrants and surfactants that bond to organic matter. This process encapsulates the bioburden so that it does not resurface on other instruments, enabling instruments to rinse freely without a greasy residue. With EcoZyme, you won't have to worry about skin irritation, respiratory issues, or that awful smell that is common among other brands.

EcoZyme Ultra Designed for: Manual Cleaning Ultrasonics Automated Washers



BEFORE



AFTER

PART ID: #7165068-SM

DoseTech Detergent Dosing Pump

1 Unit Each

*Shown in Picture

†Available for consignment or purchase



www.southmedic.com