


CO₂ Sampling and Oxygen Delivery Cannula



- Side Stream CO₂ Sampling with simultaneous oxygen delivery
- Soft anatomically curved nasal prongs for a comfortable fit
- Anatomically designed block provides a stable fit without twisting
- CO₂ is collected at the nares with a supplemental port to collect orally, while supplemental O₂ is delivered upwards to the nostrils
- 25 per box

Code	Description	Size	Tubing Length	Luer Type
NC-2025-7	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and CO ₂ connector female luer	Adult	7'	Female
NC-2025-8F	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line female luer	Adult	7'	Female
NC-2025-8	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line male luer	Adult	7'	Male
NCD-2025-8	CO ₂ /O ₂ divided cannula with 7' O ₂ with CO ₂ lines	Adult	7'	Male
NC-2025-14	CO ₂ /O ₂ nasal cannula with 14' O ₂ with universal connector and 14' CO ₂ sampling line male luer	Adult	14'	Male
NCP-2025-8	CO ₂ /O ₂ nasal cannula with 7' O ₂ with universal connector and 7' CO ₂ sampling line male luer	Paediatric	7'	Male



Respiratory

Surgical

Blades & Scalpels

Dynamic Tissue Systems

Manufactured by:
Southmedic Inc.
50 Alliance Blvd.,
Barrie, ON
Canada L4M 5K3
ISO 13485

www.southmedic.com
email: contactus@southmedic.com
1-800-463-7146 Toll free in North America
1-705-726-9383
1-705-728-9537 Fax