

# OxyArm™

- Unique oxygen delivery device with minimal contact and no facial contact
- Ideal for burn units, trauma centers and homecare
- Delivers 24% to 50% oxygen on flows from 1 to 12-15 litres per minute
- Designed for both nose and mouth breathers
- Adapts for left or right sided positioning
- Not made with natural rubber latex
- 25/case (unless otherwise stated)



Code	Description	Design	Tubing Length
OA-DA-1025	OxyArm™ Disposable Diffuser Arm only	Arm only	none
OA-DA-3025	OxyArm™ Disposable Cannula Arm only	Arm only	none
OA-PLUS-1125-8	OxyArm™ Plus Clear Adjustable Head Band with Diffuser Arm with 7' tubing	Adjustable	7'
OA-PLUS-3125-8	OxyArm™ Plus Clear Adjustable Head Band with Nasal Cannula Arm & 7" tubing	Adjustable	7"
OA-HB-1125	OxyArm™ Plus Clear Adjustable Head Band only	Adjustable	none

**Southmedic** 

- Respiratory
- Surgical
- Blades & Scalpels
- Dynamic Tissue Systems

Manufactured by:  
 Southmedic Inc.  
 50 Alliance Blvd.,  
 Barrie, ON Canada L4M 5K3  
[www.thebetteroxygenmask.com](http://www.thebetteroxygenmask.com)

email: [oxymask@southmedic.com](mailto:oxymask@southmedic.com)  
 1-800-463-7146 Toll free in North America  
 1-705-726-9383  
 1-705-728-9537 Fax  
 ISO 13485