


# OxyMask™ Aerosol




- Right angle swivel adapter allows for optimal aerosol delivery when patient is laying down or sitting
- Tapered design allows interfacing with a wide variety of nebulizer cups
- Provides directional flow to mouth and nose
- Closed sides and open bottom design:
  - prevents aerosol from going into the eyes
  - minimizes environmental contamination
  - minimizes exposure to caregiver and other patients
- Four sizes available:
  - Adult
  - Kid (3 - 10 years old, 33 - 70 lbs / 15 - 32 kg)
  - Tyke (6 months - 3 years old, 16 - 33 lbs / 7 - 15 kg)
- Not made with natural rubber latex or phthalates such as DEHP
- 25 mask/case



Aerosol flow WITHOUT obstruction:



Competitors masks express aerosol as shown upward to the eyes  
Note: tested at 8 LPM



The OxyMask Aerosol Series express aerosol towards patients nose and mouth

Aerosol flow WITH obstruction:




Competitors masks express aerosol as shown upward to the eyes



The OxyMask Aerosol Series express aerosol downward

Code	Description	Mask Size	Tubing Length	Elastic
OMN-1025	OxyMask™ Aerosol with 22 mm elbow	Adult	None	Single
OKN-1025	OxyMask™ Aerosol Kid with 22 mm elbow	Kid	None	Single
OTN-1025	OxyMask™ Aerosol Tyke with 22 mm elbow	Tyke	None	Single



Manufactured by:  
Southmedic Inc.  
50 Alliance Blvd.,  
Barrie, ON Canada L4M 5K3  
www.thebetteroxygenmask.com

email: [oxymask@southmedic.com](mailto:oxymask@southmedic.com)  
1-800-463-7146 Toll free in North America  
1-705-726-9383  
1-705-728-9537 Fax  
ISO 13485